

Supplementary Figure S1. anti-HA immunoblot for the SH-EGFR del TAP (half pulldown series). (A) 100 μ g whole cell lysate; (B) 20 mg input 1-step biotin elution; (C) 20 mg input 2-step formic acid elution; (D) 12.5 mg input 1-step biotin elution; (E) 12.5 mg input 2-step formic acid elution; (F) 5 mg input 1-step biotin elution; (G) 5 mg input 2-step formic acid elution; (H) 2.5 mg input 1-step biotin elution; and (I) 2.5 mg input 2-step formic acid elution. 6% of the biotin and final eluate were loaded. The red arrow indicates the position of tagged-SH-EGFR del.

